

O.M Shade Guide - ZiOx opak

Vita Classic	3D Master	Chromascope	Biodent	Myerson	O.M Opak
A1	1M2	01/110	13	58	CT1
A2	2M2	2A/130	20	65	CT1/CT2
A3	3L1.53M3	5B/320	23	65	CT5
A3,5	3M3	1E/230	30	65	CT7
A4	4L2.5	3C/530	32	C	CT9
B1	2M1	01/110	10	58	CT1
B2	2L1.5	2A/130	17	H	CT5
B3	3L2.5	3A/310	26	E	CT7
B4	3L2.5	5B/320	26	41	CT7
C1	2L1.5	1A/120	13	H	CT5
C2	3L1.5	4B/430	40	G/H	CT7
C3	4L1.5	6D/510	40	69	CT9
C4	5M1	4C/520	41	41	CT9
D2	3M1	4A/410	21	60	CT5
D3	3R1.5	6B/420	16	G	CT5/CT7
D4	3L1.5	6C440	30	F	CT7

O.M Shade Guide - ZiOx translucent

Vita Classic	3D Master	Chromascope	Biodent	Myerson	O.M HT	O.M Shine
A1	1M2	01/110	13	58	LL1	
A2	2M2	2A/130	20	65	LL2	X
A3	3L1.53M3	5B/320	23	65	LL2	X
A3,5	3M3	1E/230	30	65	LL3/LE2	X
A4	4L2.5	3C/530	32	C	LL4/LE3	
B1	2M1	01/110	10	58	LL1	
B2	2L1.5	2A/130	17	H	LL1	X
B3	3L2.5	3A/310	26	E	LE2	X
B4	3L2.5	5B/320	26	41	LL2/LE2	
C1	2L1.5	1A/120	13	H	LL2	
C2	3L1.5	4B/430	40	G/H	LE2	X
C3	4L1.5	6D/510	40	69	LL2/LL3	
C4	5M1	4C/520	41	41	LL4	
D2	3M1	4A/410	21	60	LL2	
D3	3R1.5	6B/420	16	G	LE2	
D4	3L1.5	6C440	30	F	LE3	

OBS! Översättningen är ungefärlig och kan endast ses som rådgivande.