

Only for 3shape user Designation	Produced by OpenMill*			Material		Scanbody									Dim analog*3
	Size	I-Bridge X	I-Butment X	CoCr	Titan	Original	Kulzer SCAB	NT-Trading*2	DAS IO	DAS Lab	Elos IO	Elos LAB	Medentika	3shape	
Ankylos Balance Base C Abut. Level	4,2	-	-	x	x	-	-	-	-	-	-	x	-	-	-
Astra Uni Abutment Level	20°	3 days	-	x	x	-	x	-	-	-	x	-	-	-	DAS
	45°	-	-	x	x	-	x	-	-	-	x	-	-	-	-
Astra implant level EV	EV	3 days	-	x	x	-	x	-	x	x	x	-	-	-	DAS
	3,0	-	-	-	x	-	-	x	x	x	x	-	x	x	DAS
	3,6	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
	4,2	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
	4,8	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
Bego Semados Implant Level	5,4	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
	3,75	-	-	x	x	x	x	x	-	-	-	-	-	-	NT
	4,1	-	-	x	x	x	x	x	-	-	-	-	x	-	NT
BioHorizons Internal Implant Level	4,5	-	-	x	x	x	x	x	-	-	-	-	x	-	NT
	5,5	-	-	x	x	x	x	x	-	-	-	-	x	-	NT
	3,0	-	-	x	x	-	-	-	-	-	-	-	-	-	-
Biomet 3i Certain Implant Level	3,5	-	-	x	x	-	x	x	-	-	x	x	-	x	NT
	4,5	-	-	x	x	-	x	x	-	-	x	x	-	x	NT
	5,7	-	-	x	x	-	x	x	-	-	x	x	-	x	NT
Biomet 3i External Hex Implant Level	3,4	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	4,1	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	5,0	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
Biomet 3i Low Profile Abutment Level	3,4	-	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	4,1	-	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	5,0	-	-	x	x	-	x	x	x	x	x	x	x	x	DAS
Nobel Biocare Brånemark Implant Level	6,0	-	-	x	x	-	x	x	-	-	x	x	-	x	-
	LP	-	-	x	x	-	x	-	x	x	x	x	-	-	DAS
	NP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
Nobel Biocare Multi Unit Abut. Level	RP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	WP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	NP/RP	3 days	-	x	x	-	x	-	x	-	x	x	-	-	DAS
Altatec Camlog Implant Level	WP	-	-	x	x	-	x	-	-	-	-	-	-	-	-
	3,3	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	3,8	-	-	x	x	-	x	x	x	x	x	x	x	x	NT
	4,3	-	-	x	x	-	x	x	x	x	x	x	x	x	NT
	5,0	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
Altatec Conelog Implant Level	6,0	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	3,3	-	-	x	x	-	x	x	x	x	-	-	x	-	NT
	3,8	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
Dentaurum tiologic Implant Level	4,3	-	-	x	x	-	x	x	x	x	-	-	x	-	NT
	5,0	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	S	-	-	x	x	-	-	-	-	-	-	-	-	-	-
Dentaurum tiologic Abutment Level	M	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	L	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	S/M/L	-	-	x	x	-	-	-	-	-	-	-	-	-	-
Medentis ICX-templant Implant Level	3,3 - 4,8	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
Neoss ProActive Implant Level	3,5 - 5,5	-	-	x	x	-	x	x	-	-	x	x	-	-	NT
Nobel Biocare NobelReplace Implant Level	NP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	RP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	WP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
Nobel Biocare NobelActive Implant Level	6,0	3 days	-	x	x	-	-	-	-	-	-	-	-	-	NT
	3,5	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	4,3	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
Straumann Bone Level Implant Level	NNC	-	-	-	-	x	-	-	x	x	x	x	x	x	DAS
	NC 3,3	-	-	-	-	x	x	x	x	x	x	x	x	x	DAS
	RC 4,1	-	-	-	-	x	x	x	x	x	x	x	x	x	DAS
Straumann Bone Level Abutment Level	RC 4,8	-	-	-	-	x	x	x	x	x	x	x	x	x	DAS
	3,5	-	-	x	x	x	-	-	-	-	-	-	-	-	-
	4,6	-	-	x	x	x	-	-	x	x	x	x	-	x	DAS
Straumann Tissue Level Implant Level	NN	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	RN	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
	WN	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
Straumann synOcta Abutment Level	RN	-	3 days	x	x	-	x	-	x	x	-	-	-	-	DAS
	WN	-	-	x	x	-	x	-	-	-	-	-	-	-	-
Zimmer Tapered Screw-Vent Implant Level	3,5	-	-	x	x	-	x	x	x	x	x	x	x	x	NT
	4,5	-	-	x	x	-	x	x	x	x	x	x	x	x	NT
	5,7	-	-	x	x	-	x	x	x	x	x	x	x	x	NT
Friadent Xive Implant Level	3,0	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	3,4	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	3,8	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	4,5	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	5,5	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
Thommen SPI Implant Level	3,5	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	4,0	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	4,5	-	-	x	x	-	-	x	-	-	-	-	-	-	NT
	5,0	-	-	x	x	-	-	x	-	-	-	-	-	-	NT
Thommen VARIOmulti Abutment Level	6,0	-	-	x	x	-	-	x	-	-	-	-	-	-	NT
	-	-	-	x	x	-	x	-	-	-	-	-	-	-	NT
Sweden & Martina Premium Kohno Implant Level	3,3	x	x	x	x	-	x	x	x	x	-	-	-	-	NT
	3,8	x	x	x	x	-	x	x	-	-	-	-	-	-	NT
	4,25	x	x	x	x	-	x	x	-	-	-	-	-	-	NT
	5	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
Bredent SKY Implant Level	3,5 - 4,0	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
OSstem Multi Unit Abutment	NP/RP	3 days	-	x	x	-	x	-	x	-	x	x	-	-	DAS

Produced by OpenMill* = 3 days means we produce it in-house in 3 days.

NT Trading*2 = Scanbodies and other components from NT-

Dim analog*3 = DIM stands for Digital Implant Model and are analogs used primarily for digitally 3d-printed models.