

Designation	Produced by OpenMill*		Material		Scanbody										Dim analog*3
	Size	I-Bridge X	I-Butment X	CoCr	Titan	Original	Kulzer SCAB	NT-Trading*2	DAS IO	DAS Lab	Elos IO	Elos LAB	Medentika	3shape	
Ankylos Balance Base C Abut. Level	4.2	-	-	x	x	-	-	-	-	-	-	x	-	-	-
Astra Uni Abutment Level	20°	3 days	-	x	x	-	x	-	-	-	x	-	-	-	-
	45°	-	-	x	x	-	x	-	-	-	x	-	-	-	-
	EV	3 days	-	x	x	-	x	-	x	x	x	-	-	-	DAS
Astra implant level EV	3,6	-	-	x	x	-	-	x	x	x	x	-	x	x	DAS
	4,2	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
	4,8	-	-	x	x	-	-	x	x	x	x	-	x	x	DAS
	5,4	-	-	x	x	-	-	x	x	x	x	-	x	x	DAS
Bego Semados Implant Level	3.75	-	-	x	x	x	x	x	-	-	-	-	-	-	NT
	4.1	-	-	x	x	x	x	x	-	-	-	-	x	-	NT
	4.5	-	-	x	x	x	x	x	-	-	-	-	x	-	NT
	5.5	-	-	x	x	x	x	x	-	-	-	-	x	-	NT
BioHorizons Internal Implant Level	3.0	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	3.5	-	-	x	x	-	x	x	-	-	x	x	-	x	NT
	4.5	-	-	x	x	-	x	x	-	-	x	x	-	x	NT
	5.7	-	-	x	x	-	x	x	-	-	x	x	-	x	NT
Biomet 3i Certain Implant Level	3,4	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	4,1	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	5,0	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
Biomet 3i External Hex Implant Level	3,4	-	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	4,1	-	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	5,0	-	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	6.0	-	-	x	x	-	x	x	-	-	x	x	-	x	-
Biomet 3i Low Profile Abutment Level	-	-	-	x	x	-	x	-	-	-	-	-	-	-	-
Nobel Biocare Brånemark Implant Level	NP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	RP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	WP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
Nobel Biocare Multi Unit Abut. Level	NP/RP	3 days	-	x	x	-	x	-	x	-	x	x	-	-	DAS
	WP	-	-	x	x	-	x	-	-	-	-	-	-	-	-
Altatec Camlog Implant Level	3.3	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	3.8	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	4.3	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	5.0	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	6.0	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
Altatec Conelog Implant Level	3.3	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
	3.8	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
	4.3	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
	5.0	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
Dentaurum tiologic Implant Level	S	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	M	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	L	-	-	x	x	-	-	-	-	-	-	-	-	-	-
Dentaurum tiologic Abutment Level	S/M/L	-	-	x	x	-	-	-	-	-	-	-	-	-	-
Medentis ICX-templant Implant Level	3.3 - 4.8	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
Neoss ProActive Implant Level	3.5 - 5.5	-	-	x	x	-	x	x	-	-	x	x	-	-	NT
Nobel Biocare NobelReplace Implant Level	NP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	RP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	WP	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	6,0	3 days	-	x	x	-	-	-	-	-	-	-	-	-	-
Nobel Biocare NobelActive Implant Level	3,5	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	4.3	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
Straumann Bone Level Implant Level	NNC	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	NC 3,3	-	-	-	-	-	x	x	x	x	x	x	x	x	DAS
	RC 4.1	-	-	-	-	-	x	x	x	x	x	x	x	x	DAS
	RC 4.8	-	-	-	-	-	-	-	-	-	-	-	x	-	-
Straumann Bone Level Abutment Level	3.5	-	-	x	x	x	-	-	-	-	-	-	-	-	-
	4.6	-	-	x	x	x	-	-	x	x	x	x	-	x	DAS
Straumann Tissue Level Implant Level	NN	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	RN	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
	WN	3 days	3 days	x	x	-	x	x	x	-	x	x	x	x	DAS
Straumann synOcta Abutment Level	RN	-	-	x	x	-	x	-	x	x	-	-	-	-	DAS
	WN	-	-	x	x	-	x	-	-	-	-	-	-	-	-
Zimmer Tapered Screw-Vent Implant Level	3,5	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	4,5	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	5,7	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
Friadent Xive Implant Level	3,0	-	-	x	x	-	-	-	-	-	-	-	-	-	-
	3,4	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	3,8	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	4,5	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
	5,5	-	-	x	x	-	x	x	-	-	x	x	x	x	NT
Thommen SPI Implant Level	3,5	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	4,0	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	4,5	-	-	x	x	-	-	x	-	-	-	-	-	-	NT
	5,0	-	-	x	x	-	-	x	-	-	-	-	-	-	NT
	6,0	-	-	x	x	-	-	x	-	-	-	-	-	-	NT
Thommen VARIOmulti Abutment Level	-	-	-	x	x	-	x	-	-	-	-	-	-	-	NT
Sweden & Martina Premium Kohno Implant Level	3,3	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	3,8	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	4,25	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
	5	-	-	x	x	-	x	x	-	-	-	-	-	-	NT
Bredent SKY Implant Level	3.5 - 4,0	-	-	x	x	-	x	x	-	-	-	-	x	-	NT
OSstem Multi Unit Abutment	NP/RP	3 days	-	x	x	-	x	-	x	-	x	x	-	-	-

Produced by OpenMill* = 3 days means we produce it in-house in 3 days.

NT Trading*2 = Scanbodies and other components from

Dim analog*3 = DIM stands for Digital Implant Model and are analogs used primarily for digitally 3d-printed models.