

Only for 3shape user	Produced by OpenMill*			Material		Scanbody									Dim analog* ³
Designation	Size	I-Bridge X	I-Butment X	CoCr	Titan	Original	Kulzer SCAB	NT-Trading* ²	DAS IO	DAS Lab	Elos IO	Elos LAB	Medentika	3shape	
Astra Uni Abutment Level	20°	3 days	-	x	x	-	x	-	-	-	x	-	-	-	-
	45°	-	-	x	x	-	x	-	-	-	x	-	-	-	-
	EV	3 days	-	x	x	-	x	-	x	x	x	-	-	-	DAS
Astra implant level EV	3,6	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
	4,2	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
	4,8	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
	5,4	-	3 days	x	x	-	-	x	x	x	x	-	x	x	DAS
Biomet 3i Certain Implant Level	3,4	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	4,1	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
	5,0	3 days	3 days	x	x	-	x	x	x	x	x	x	x	x	DAS
Biomet 3i External Hex Implant Level	3,4	3 days	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	4,1	3 days	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	5.0	3 days	-	x	x	-	x	x	x	x	x	x	x	x	DAS
	6.0	3 days	-	x	x	-	x	x	-	-	x	x	-	x	-
Nobel Biocare Brånemark Implant Level	NP	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
	RP	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
	WP	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
Nobel Biocare Multi Unit Abut. Level	NP/RP	3 days	-	x	x	x	x	-	x	-	x	x	-	-	DAS
	WP	NA	-	x	x	x	x	-	-	-	-	-	-	-	-
Neoss ProActive Implant Level	3.5 - 5.5	3 days	NA	x	x	-	x	x	-	-	x	x	-	-	NT
Nobel Biocare NobelReplace Implant Level	NP	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
	RP	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
	WP	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
	6,0	3 days	NA	x	x	x	-	-	-	-	-	-	-	-	-
Nobel Biocare NobelActive Implant Level	3,5	3 days	3 days	x	x	x	x	x	x	x	x	x	x	x	DAS
	4.3	3 days	3 days	x	x	x	x	x	x	x	x	x	x	x	DAS
Straumann Bone Level Implant Level	NC	NA	-	-	-	-	-	-	-	-	-	-	-	-	-
	NC 3,3	3 days	-	-	-	x	x	x	x	x	x	x	x	x	DAS
	RC 4.1/4.8	3 days	-	-	-	x	x	x	x	x	x	x	x	x	DAS
Straumann Bone Level Abutment Level	3.5	3 days	NA	x	x	x	-	-	-	-	-	-	-	-	-
	4.6	3 days	NA	x	x	x	-	-	x	x	x	x	-	x	DAS
Straumann Tissue Level Implant Level	NNC	NA	NA	x	x	-	-	-	-	-	-	-	-	-	-
	RN	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
	WN	3 days	3 days	x	x	x	x	x	x	-	x	x	x	x	DAS
Straumann synOcta Abutment Level	RN	3 days	3 days	x	x	-	x	-	x	x	-	-	-	-	DAS
	WN	3 days	3 days	x	x	-	x	-	-	-	-	-	-	-	-
OSstem Multi Unit Abutment	NP/RP	3 days	NA	x	x	-	x	-	x	-	x	x	-	-	-
Produced by OpenMill* = 3 days means we produce it in-house in 3 days.		NT Trading* ² = Scanbodies and other components from			Dim analog* ³ = DIM stands for Digital Implant Model and are analogs used primarily for digitally 3d-printed models.										